

momo therapy bicycle.

Firmly in the saddle.

GMFCS
Level I-V

kg
max. 100



30-95 cm

16-23 cm

I am your therapy bicycle for more **balance.**

I am there for you if the pleasant weather tempts you to take a trip around the lake. But I can do much more than that if you take me on an excursion. Depending on how steady you are on your bike, my stabilizer wheels can be adjusted from rigid to very soft spring settings. In this way, you can easily train your balance. Apart from that, I look really cool: without my supporting functions, I look just like a bicycle!





momo therapy bicycle.

Handlebars

Available in many different designs.

Individual drive

from a rigid drive via a 7-gear freewheel brake hub, right up to electrical support.

Saddle

available in different shapes and sizes.

Two brake systems

from 16", consisting of rims and drum brakes.

Pedals

can be individually configured.

Stabilizer wheels

can be adjusted from rigid to very soft.

Low access level as standard feature

permits the rider to climb on conveniently and also allows independent transfer.

You can find a diverse range of accessories for optimum support & individual equipment from page 132.

You have the choice! Between many handlebars, saddles, foot pans, head supports and different drives - from rigid toothed wheels right up to the 3-gear hub gear with reverse option!



Frame colours 12"

- Telemagenta
- Sky blue
- Pure white

Frame colours 16" and 20"

- Melon yellow
- Navy blue
- Pure white

Frame colours 24" and 26"

- Anthracite metallic
- Pure white




Using me, you can train your balance, coordination and endurance.



Dynamic stabilizer wheels

can be adjusted to your capabilities. If set to rigid, they will initially give you security and build trust. Over the course of your therapy, they can be adjusted to become very gentle suspension stabilizers. These can easily balance you as you go round bends and also offer sufficient drive force on dirt roads.



Data	12"	16"	20"	24"	26"
Pedal to top edge of saddle*	30 - 38 cm	42 - 57,5 cm	50 - 64,5 cm	59 - 76,5 cm	71,5 - 87,5 cm
Pedal to top edge of saddle**	33 - 41 cm	45 - 60 cm	52 - 68 cm	61 - 84 cm	73 - 95 cm
Saddle support tube to front tube	33 cm	39 cm	45 cm	53 cm	58 cm
Wheelbase	64 cm	76 cm	89 cm	104 cm	112 cm
Full length	99 cm	122 cm	148 cm	168 cm	178 cm
Access point height	16 cm	18 cm	20 cm	22 cm	23 cm
Full width	43 cm	52 cm	62 cm	71 cm	76 cm
max. load	35 kg	50 kg	65 kg	100 kg	100 kg
Weight	10,5 kg	13 kg	15 kg	18 kg	20 kg
HMV No.	22.51.01.0014				

* Dimensions with standard saddle post / ** Dimensions with T-saddle post



You can find further information on accessories & harnesses
on pages 132 - 135 and in the price list.



momo therapy bicycle.

With E-drive.

Did you know that your momo therapy bicycle. can also be equipped with an electric drive? This supports you where your own power is not quite sufficient. Whether it's a short or long trip to town or to the lake.

The basis of our system is Heinzmann's DirectPower e-drive in the front wheel. This is a noiseless direct drive, which is optionally available with a starting aid. For experienced drivers the drive is also available in an "advanced version". This allows a more powerful drive. The "beginner's version" can also be used to control the Eco, Standard or Power mode via display.

With this power-supporting electric drive we offer an ideal addition to your **momo therapy bicycle.** and help you to achieve your set targets.

We are at your disposal for an individual consultation.

DirectPower E-Bike drive

is designed brushless and without gears, minimizes maintenance and wear and tear and eliminates motor noise.



Button or twist handle operation

allow you to set the mode, and the start-up assistance if necessary, including during travel.



LCD display

shows you information on the support level, speed and battery status - and is easy to read, even in bright / sunny surroundings.





The battery can also be removed for charging. The charge level can be read on the display.

