
Marcy. There for you from the start.




I am your anterior walking aid, and will assist you in your very first attempts


to walk. I can hardly wait to start with you! I will support your balance and simultaneously train your torso musculature. You can adapt my wheelbase to suit: if you adjust it so that it is somewhat shorter, I become highly manoeuvrable and you can move me better. I give you more security and support with a long wheelbase. If you find it hard to walk in a straight line, I am happy to assist you with my direction locks.



















Data	Size O	Size 1	Size 2	Size 3
Handle height	30 - 47 cm	38 - 59 cm	53 - 75 cm	63 - 88 cm
Height of forearm rests	41 - 58 cm	49 - 70 cm	63 - 85cm	73 - 98 cm
Handle spacing	18 - 32 cm	22 - 36 cm	30 - 44 cm	36 - 50 cm
Seat height	25 cm	30 cm	45 cm	45 cm
Total width	44 cm	48 cm	62 cm	71cm
Total length	36,5 cm	42,5 cm	57,5 cm	61cm
Total length (telescopic wheelbase)	39 - 47 cm	45 - 55cm	60 - 74 cm	65 - 79cm
max. load	20 kg	30 kg	50 kg	65 kg
Weight	4,2 kg	5 kg	7,4 kg	9,1 kg
HMV No.	10.46.02.3049			





