
rev. 01- 16.02.23



by Orthosanit Diffusion

Main features

Comfort

Pad is made of high quality foam, with a removable and washable cover realized with technical material for maximal comfort.

Simplicity

Smooth and simple style combined with black color merge into beauty and ergonomics.

Strength

Being a single piece of aluminium without any other separated parts it guarantees a long and lasting safety

Lightness

Thanks to minimum aluminium thickness, Active allows to keep the wheelchair's wheight as low as possible



Size



 LOW
 33 cm

 MID
 41 cm

 HIGH
 48 cm



Swing-Away Lateral support

- **S** for wheelchairs from 30 to 35
- $\ \ \, \mathbf{M}$ for wheelchairs from 36 to 40
- **L** for wheelchairs from 41 to 45
- XL for wheelchairs from 46 to 49

By combining and adjusting the thoracal supports you can provide the right control and support tu the user





ACTIVE – Height/ Size /Lateral support Page 2











3 heights 4 sizes


10 models




By adding the lateral supports you can reach maximum customization

Adjustable clamp for rigid frame wheelchairs









ACTIVE - Adjustable clamp for rigid frame wheelchairs Page 4


Adjustable clamp for folding wheelchairs












Features:

- 1 rigidizer bar to specify:
 - **S**: 30 36 cm / **M**: 36 41 cm / **L**: 41 46 cm
- 2 hose clamp: 19-21 / 22-25 mm
- (3) height adjustment
- depth adjustment 4
- angle adjustment 5










ACTIVE - Adjustable clamp for folding wheelchairs Page 5





Balance points

on the rigidizer bar to simplify the balance of the clamp for folding wheelchairs.

2 rubber tubes (internal diameter 22mm) easy adaptables to every wheelchair **to protect the backrest.**





Highly breathable pads removable and **washable** at 30°.

